

产品名称: SNX6 Mouse Monoclonal Antibody
产品货号: AMM85922



产品概述 (Summary)

产品名称 (Production Name)	SNX6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX6
别名 (Alternative Names)	Sorting nexin-6, TRAF4-associated factor 2, Sorting nexin-6, N-terminally processed, SNX6
基因 ID (Gene ID)	58533.0
蛋白 ID (SwissProt ID)	Q9UNH7. This SNX6 monoclonal antibody is generated from mouse immunized with SNX6 recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:15000
蛋白分子量 (Molecular Weight)	46.6kDa

研究背景 (Background)


产品名称: SNX6 Mouse Monoclonal Antibody
产品货号: AMM85922



This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein associates with the long isoform of the leptin receptor, the transforming growth factor-beta family of receptor serine-threonine kinases, and with receptor tyrosine kinases for platelet-derived growth factor, insulin, and epidermal growth factor. This protein may form oligomeric complexes with family member proteins through interactions of both the PX domain and the coiled coil regions of the molecules. Translocation of this protein from the cytoplasm to the nucleus occurs after binding to proviral integration site 1 protein. This gene results in two transcripts encoding two distinct isoforms.

研究领域 (Research Area)

图片 (Image Data)



SNX6 Antibody (Cat. #AMM85922) western blot analysis in U251 cell line lysates (35µg/lane). This demonstrates the SNX6 antibody detected the SNX6 protein (arrow).

注意事项 (Note)

For research use only .