

产品名称: GBX2(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85924



产品概述 (Summary)

产品名称 (Production Name)	GBX2(Ascites) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgM
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GBX2(Ascites)
别名 (Alternative Names)	Homeobox protein GBX-2, Gastrulation and brain-specific homeobox protein 2, GBX2
基因 ID (Gene ID)	2637.0
蛋白 ID (SwissProt ID)	P52951. This GBX2 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 103-131 amino acids from the Central region of human GBX2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:4000-1:64000
蛋白分子量 (Molecular Weight)	37.3kDa

产品名称: GBX2(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85924




研究背景 (Background)

GBX2 may act as a transcription factor for cell pluripotency and differentiation in the embryo.

研究领域 (Research Area)

图片 (Image Data)



GBX2 Antibody (Center) (Cat. #AMM85924) western blot analysis in human placenta tissue lysates (35 μ g/lane). This demonstrates the GBX2 antibody detected the GBX2 protein (arrow).

注意事项 (Note)

For research use only .