

产品名称: NCOA4(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85928



产品概述 (Summary)

产品名称 (Production Name)	NCOA4(Ascites) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgM
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCOA4(Ascites)
别名 (Alternative Names)	Nuclear receptor coactivator 4, NCoA-4, Androgen receptor coactivator 70 kDa protein, 70 kDa AR-activator, 70 kDa androgen receptor coactivator, Androgen receptor-associated protein of 70 kDa, Ret-activating protein ELE1, NCOA4, ARA70, ELE1, RFG
基因 ID (Gene ID)	8031.0
蛋白 ID (SwissProt ID)	Q13772.Purified His-tagged NCOA4 protein(Fragment) was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:3200
蛋白分子量 (Molecular Weight)	69.7kDa

产品名称: NCOA4(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85928




研究背景 (Background)

This gene encodes an androgen receptor coactivator. The encoded protein interacts with the androgen receptor in a ligand-dependent manner to enhance its transcriptional activity. Chromosomal translocations between this gene and the ret tyrosine kinase gene, also located on chromosome 10, have been associated with papillary thyroid carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes are present on chromosomes 4, 5, 10, and 14.

研究领域 (Research Area)

图片 (Image Data)



NCOA4 Antibody (Cat. #AMM85928) western blot analysis in 293 cell line lysates (35µg/lane). This demonstrates the NCOA4 antibody detected the NCOA4 protein (arrow).

注意事项 (Note)

For research use only .