

产品名称: DYT1 Mouse Monoclonal Antibody
产品货号: AMM85937



产品概述 (Summary)

产品名称 (Production Name)	DYT1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DYT1
别名 (Alternative Names)	Torsin-1A, Dystonia 1 protein, Torsin ATPase-1A, 364-, Torsin family 1 member A, TOR1A, DQ2, DYT1, TA, TORA
基因 ID (Gene ID)	1861.0
蛋白 ID (SwissProt ID)	O14656. Purified His-tagged DYT1 protein(Fragment) was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:8000
蛋白分子量 (Molecular Weight)	37.9kDa


研究背景 (Background)

产品名称: DYT1 Mouse Monoclonal Antibody
产品货号: AMM85937

The protein encoded by this gene is a member of the AAA family of adenosine triphosphatases (ATPases), is related to the Clp protease/heat shock family and is expressed prominently in the substantia nigra pars compacta. Mutations in this gene result in the autosomal dominant disorder, torsion dystonia 1. [provided by RefSeq].

研究领域 (Research Area)

图片 (Image Data)



DYT1 Antibody (Cat. #AMM85937) western blot analysis in ZR-75-1 cell line lysates (35µg/lane). This demonstrates the DYT1 antibody detected the DYT1 protein (arrow).

注意事项 (Note)

For research use only .