

产品名称: CD3 Mouse Monoclonal Antibody
产品货号: AMM85941



产品概述 (Summary)

产品名称 (Production Name)	CD3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD3
别名 (Alternative Names)	T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T3E
基因 ID (Gene ID)	916.0
蛋白 ID (SwissProt ID)	P07766. This CD3 Monoclonal antibody is generated from mice immunized with a KLH conjugated synthetic peptide selected from human CD3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:200-1:1600
蛋白分子量 (Molecular Weight)	23.1kDa

研究背景 (Background)


产品名称: CD3 Mouse Monoclonal Antibody
产品货号: AMM85941



The CD3 complex mediates signal transduction.

研究领域 (Research Area)

图片 (Image Data)



CD3 Antibody (Ascites)(Cat. #AMM85941) western blot analysis in mouse spleen tissue lysates (35 μ g/lane). This demonstrates the CD3 antibody detected the CD3 protein (arrow).

注意事项 (Note)

For research use only .