

产品名称: FN3K (N-term)(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85942



产品概述 (Summary)

产品名称 (Production Name)	FN3K (N-term)(Ascites) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FN3K (N-term)(Ascites)
别名 (Alternative Names)	Fructosamine-3-kinase, 271-, FN3K
基因 ID (Gene ID)	64122.0
蛋白 ID (SwissProt ID)	Q9H479. Purified His-tagged FN3K protein was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:8000
蛋白分子量 (Molecular Weight)	35.2kDa

研究背景 (Background)


产品名称: FN3K (N-term)(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85942



May initiate a process leading to the deglycation of fructoselysine and of glycated proteins. May play a role in the phosphorylation of 1-deoxy-1-morpholinofructose (DMF), fructoselysine, fructoseglycine, fructose and glycated lysozyme.

研究领域 (Research Area)

图片 (Image Data)



FN3K Antibody (N-term) (Cat. #AMM85942) western blot analysis in K562 cell line lysates (35µg/lane). This demonstrates the FN3K antibody detected the FN3K protein (arrow).

注意事项 (Note)

For research use only .