

产品名称: CCRK (N-term) Mouse Monoclonal Antibody
产品货号: AMM85944



产品概述 (Summary)

产品名称 (Production Name)	CCRK (N-term) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCRK (N-term)
别名 (Alternative Names)	Cyclin-dependent kinase 20, CDK-activating kinase p42, CAK-kinase p42, Cell cycle-related kinase, Cell division protein kinase 20, Cyclin-dependent protein kinase H, Cyclin-kinase-activating kinase p42, CDK20, CCRK, CDCH
基因 ID (Gene ID)	23552.0
蛋白 ID (SwissProt ID)	Q8IZL9.Purified His-tagged CCRK protein was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	38.7kDa

产品名称: CCRK (N-term) Mouse Monoclonal Antibody
产品货号: AMM85944




研究背景 (Background)

Required for high-level Shh responses in the developing neural tube. Together with BROMI, controls the structure of the primary cilium by coordinating assembly of the ciliary membrane and axoneme, allowing GLI2 to be properly activated in response to SHH signaling (By similarity). Involved in cell growth. Activates CDK2, a kinase involved in the control of the cell cycle, by phosphorylating residue 'Thr-160'.

研究领域 (Research Area)

图片 (Image Data)



CCRK Antibody (N-term)(Cat. #AMM85944) western blot analysis in MCF-7,A549 cell line lysates (35µg/lane).This demonstrates the CCRK antibody detected the CCRK protein (arrow).

注意事项 (Note)

For research use only .