

产品名称: FAPB3(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85945



产品概述 (Summary)

产品名称 (Production Name)	FAPB3(Ascites) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAPB3(Ascites)
别名 (Alternative Names)	Fatty acid-binding protein, heart, Fatty acid-binding protein 3, Heart-type fatty acid-binding protein, H-FABP, Mammary-derived growth inhibitor, MDGI, Muscle fatty acid-binding protein, M-FABP, FABP3, FABP11, MDGI
基因 ID (Gene ID)	2170.0
蛋白 ID (SwissProt ID)	P05413. Purified His-tagged FAPB3 protein was used to produce this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:5000
蛋白分子量 (Molecular Weight)	14.9kDa

产品名称: FAPB3(Ascites) Mouse Monoclonal Antibody
产品货号: AMM85945




研究背景 (Background)

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

研究领域 (Research Area)

图片 (Image Data)



FAPB3 Antibody(Cat. #AMM85945) western blot analysis in mouse heart tissue lysates (35 μ g/lane).This demonstrates the FAPB3 antibody detected the FAPB3 protein (arrow).

注意事项 (Note)

For research use only .