

产品名称: FUBP3 (N-term) Mouse Monoclonal Antibody
产品货号: AMM85951



产品概述 (Summary)

产品名称 (Production Name)	FUBP3 (N-term) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUBP3 (N-term)
别名 (Alternative Names)	Far upstream element-binding protein 3, FUSE-binding protein 3, FUBP3, FBP3
基因 ID (Gene ID)	8939.0
蛋白 ID (SwissProt ID)	Q96I24. Purified His-tagged FUBP3 protein was used to produce this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC 1:25-1:50
蛋白分子量 (Molecular Weight)	61.4kDa

研究背景 (Background)


May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.

产品名称: FUBP3 (N-term) Mouse Monoclonal Antibody
产品货号: AMM85951




研究领域 (Research Area)

图片 (Image Data)



FUBP3 Antibody (N-term) (Cat. #AMM85951) western blot analysis in HeLa,U251 cell line lysates (35µg/lane).This demonstrates the FUBP3 antibody detected the FUBP3 protein (arrow).



Fluorescent image of HeLa cells stained with FUBP3 Antibody (N-term)(Cat#AMM85951). AMM85951 was diluted at 1:25 dilution. An Alexa Fluor® 488-conjugated goat anti-mouse IgG at 1:400 dilution was used as the secondary antibody (green). Cytoplasmic actin was counterstained with Alexa Fluor® 555 conjugated with Phalloidin (red).

注意事项 (Note)

For research use only .