

产品名称: CREB3L4 (M01) Mouse Monoclonal Antibody
产品货号: AMM85952



产品概述 (Summary)

产品名称 (Production Name)	CREB3L4 (M01) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CREB3L4 (M01) Cyclic AMP-responsive element-binding protein 3-like protein 4, cAMP-responsive element-binding protein 3-like protein 4, Androgen-induced basic leucine zipper protein, AlbZIP, Attaching to CRE-like 1, ATCE1, Cyclic AMP-responsive element-binding protein 4, CREB-4, cAMP-responsive element-binding protein 4, Transcript induced in spermiogenesis protein 40, Tisp40, hJAL, Processed cyclic AMP-responsive element-binding protein 3-like protein 4, CREB3L4, AIBZIP, CREB4, JAL
别名 (Alternative Names)	
基因 ID (Gene ID)	148327.0
蛋白 ID (SwissProt ID)	Q8TEY5. This CREB3L4 antibody is generated from mice immunized with a recombinant between 1-300 amino acids from the region of human CREB3L4.

产品应用 (Application)

产品名称: CREB3L4 (M01) Mouse Monoclonal Antibody
产品货号: AMM85952



稀释比 (Dilution Ratio) WB 1:500-1:2000
蛋白分子量 (Molecular Weight) 43.4kDa


研究背景 (Background)

Transcriptional activator that may play a role in the unfolded protein response. Binds to the UPR element (UPRE) but not to CRE element. Preferentially binds DNA with to the consensus sequence 5'-T[GT]ACGT[GA][GT]-3' and has transcriptional activation activity from UPRE. Binds to NF-kappa-B site and has transcriptional activation activity from NF-kappa-B-containing regulatory elements (By similarity).

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



CREB3L4 Antibody (M01) (Cat.# AMM85952) western blot analysis in Daudi,PC-3,MCF-7, rat C6 and L6 cell line ,rat brain tissue lysates (35µg/lane).This demonstrates the CREB3L4 (M01)antibody detected the CREB3L4 (M01)protein (arrow).

注意事项 (Note)

For research use only .