

产品名称: USP25 Mouse Monoclonal Antibody
产品货号: AMM85954



产品概述 (Summary)

产品名称 (Production Name)	USP25 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP25 Ubiquitin carboxyl-terminal hydrolase 25, Deubiquitinating enzyme 25, USP
别名 (Alternative Names)	on chromosome 21, Ubiquitin thioesterase 25, Ubiquitin-specific-processing protease 25, USP25, USP21
基因 ID (Gene ID)	29761.0
蛋白 ID (SwissProt ID)	Q9UHP3. This USP25 antibody is generated from a mouse immunized with a recombinant protein from human USP25.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	122.2kDa

产品名称: USP25 Mouse Monoclonal Antibody
产品货号: AMM85954




研究背景 (Background)

Deubiquitinating enzyme that hydrolyzes ubiquitin moieties conjugated to substrates and thus, functions to process newly synthesized Ubiquitin, to recycle ubiquitin molecules or to edit polyubiquitin chains and prevents proteasomal degradation of substrates. Hydrolyzes both 'Lys-48'- and 'Lys-63'-linked tetraubiquitin chains.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from Daudi cell line, using USP25 Antibody (Cat. #AMM85954). AMM85954 was diluted at 1:1000. A goat anti-mouse IgG H&L(HRP) at 1:5000 dilution was used as the secondary antibody. Lysate at 35µg.

注意事项 (Note)

For research use only .