

产品名称: Myoglobin Mouse Monoclonal Antibody
产品货号: AMM85965



产品概述 (Summary)

产品名称 (Production Name)	Myoglobin Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Myoglobin
别名 (Alternative Names)	Myoglobin, MB
基因 ID (Gene ID)	4151.0
蛋白 ID (SwissProt ID)	P02144. This Myoglobin antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 118-150 amino acids from the C-terminal region of human Myoglobin.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	17.2kDa

研究背景 (Background)


产品名称: Myoglobin Mouse Monoclonal Antibody
产品货号: AMM85965



Serves as a reserve supply of oxygen and facilitates the movement of oxygen within muscles.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from mouse heart and skeletal muscle tissue lysates (from left to right), using Myoglobin Antibody (C-term)(Cat. #AMM85965). AMM85965 was diluted at 1:1000 at each lane. A goat anti-mouse IgG H&L(HRP) at 1:3000 dilution was used as the secondary antibody. Lysates at 35 μ g per lane.

注意事项 (Note)

For research use only .