

产品名称: Ntrk3 Mouse Monoclonal Antibody
产品货号: AMM85971



产品概述 (Summary)

产品名称 (Production Name)	Ntrk3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ntrk3
别名 (Alternative Names)	NT-3 growth factor receptor, GP145-TrkC, Trk-C, Neurotrophic tyrosine kinase receptor type 3, TrkC tyrosine kinase, Ntrk3, TrkC
基因 ID (Gene ID)	18213.0
蛋白 ID (SwissProt ID)	Q6VNS1. This Ntrk3 antibody is generated from a mouse immunized with a recombinant protein of human Ntrk3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	92.8kDa

研究背景 (Background)


产品名称: Ntrk3 Mouse Monoclonal Antibody
产品货号: AMM85971




Receptor for neurotrophin-3 (NT-3). This is a tyrosine- protein kinase receptor. Known substrates for the Trk receptors are SHC1, PI 3-kinase, and PLC-gamma-1 (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from mouse brain, rat brain tissue lysate (from left to right), using Ntrk3 Antibody(Cat. #AMM85971). AMM85971 was diluted at 1:1000 at each lane. A goat anti-mouse IgG H&L(HRP) at 1:3000 dilution was used as the secondary antibody. Lysates at 20µg per lane.



Western blot analysis of lysate from human brain tissue lysate, using Ntrk3 Antibody(Cat. #AMM85971). AMM85971 was diluted at 1:1000. A goat anti-mouse IgG H&L(HRP) at 1:3000 dilution was used as the secondary antibody. Lysate at 20µg.

注意事项 (Note)

For research use only .