

产品名称: PPP2R1B Mouse Monoclonal Antibody
产品货号: AMM85987



产品概述 (Summary)

产品名称 (Production Name)	PPP2R1B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP2R1B
别名 (Alternative Names)	Serine/threonine-protein phosphatase 2A 65 kDa regulatory subunit A beta isoform, PP2A subunit A isoform PR65-beta, PP2A subunit A isoform R1-beta, PPP2R1B
基因 ID (Gene ID)	5519.0
蛋白 ID (SwissProt ID)	P30154.This PPP2R1B antibody is generated from a mouse immunized with a recombinant protein of human PPP2R1B.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:4000,IHC 1:100-1:500,ICC 1:25-1:50,FC 1:25-1:50
蛋白分子量 (Molecular Weight)	66.2kDa

产品名称: PPP2R1B Mouse Monoclonal Antibody
产品货号: AMM85987




研究背景 (Background)

The PR65 subunit of protein phosphatase 2A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit and a variable regulatory B subunit.

研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway

图片 (Image Data)



All lanes : Anti-PPP2R1B Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .