

产品名称: EEF1E1 Mouse Monoclonal Antibody
产品货号: AMM85989



产品概述 (Summary)

产品名称 (Production Name)	EEF1E1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1E1
别名 (Alternative Names)	Eukaryotic translation elongation factor 1 epsilon-1, Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3, Elongation factor p18, Multisynthase complex auxiliary component p18, EEF1E1, AIMP3, P18
基因 ID (Gene ID)	9521.0
蛋白 ID (SwissProt ID)	O43324. This EEF1E1 antibody is generated from a mouse immunized with a recombinant protein of human EEF1E1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:4000, ICC 1:25-1:50, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	19.8kDa

产品名称: EEF1E1 Mouse Monoclonal Antibody
产品货号: AMM85989




研究背景 (Background)

Positive modulator of ATM response to DNA damage.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-EEF1E1 Antibody at 1:2000-1:4000 dilution

注意事项 (Note)

For research use only .