

产品名称: DPYSL5 Mouse Monoclonal Antibody
产品货号: AMM85990



产品概述 (Summary)

产品名称 (Production Name)	DPYSL5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPYSL5
别名 (Alternative Names)	Dihydropyrimidinase-related protein 5, DRP-5, CRMP3-associated molecule, CRAM, Collapsin response mediator protein 5, CRMP-5, UNC33-like phosphoprotein 6, ULIP-6, DPYSL5, CRMP5, ULIP6
基因 ID (Gene ID)	56896.0
蛋白 ID (SwissProt ID)	Q9BPU6. This DPYSL5 antibody is generated from a mouse immunized with a recombinant protein of human DPYSL5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:25-1:50
蛋白分子量 (Molecular Weight)	61.4kDa

产品名称: DPYSL5 Mouse Monoclonal Antibody
产品货号: AMM85990




研究背景 (Background)

May have a function in neuronal differentiation and/or axon growth.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-DPYL5 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .