

产品名称: AVPR1A Mouse Monoclonal Antibody
产品货号: AMM85994



产品概述 (Summary)

产品名称 (Production Name)	AVPR1A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgM
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AVPR1A
别名 (Alternative Names)	Vasopressin V1a receptor, V1aR, AVPR V1a, Antidiuretic hormone receptor 1a, Vascular/hepatic-type arginine vasopressin receptor, AVPR1A, AVPR1
基因 ID (Gene ID)	552.0
蛋白 ID (SwissProt ID)	P37288. This AVPR1A Monoclonal antibody is generated from mice immunized with a KLH conjugated synthetic peptide selected from the 353-383 region of human AVPR1A.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	46.8kDa

产品名称: AVPR1A Mouse Monoclonal Antibody
产品货号: AMM85994




研究背景 (Background)

Receptor for arginine vasopressin. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system. Has been involved in social behaviors, including affiliation and attachment.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-AVPR1A Antibody (C-term) at 1:2000 dilution

注意事项 (Note)

For research use only .