

产品名称: ALDOC Mouse Monoclonal Antibody
产品货号: AMM85995



产品概述 (Summary)

产品名称 (Production Name)	ALDOC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG3
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALDOC
别名 (Alternative Names)	Fructose-bisphosphate aldolase C, Brain-type aldolase, ALDOC, ALDC
基因 ID (Gene ID)	230.0
蛋白 ID (SwissProt ID)	P09972. Purified His-tagged ALDOC protein was used to produce this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	39.5kDa


研究背景 (Background)

产品名称: ALDOC Mouse Monoclonal Antibody
产品货号: AMM85995



研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-ALDOC Antibody (C-term) at 1:2000 dilution

注意事项 (Note)

For research use only .