

产品名称: TUSC2 Mouse Monoclonal Antibody
产品货号: AMM86001



产品概述 (Summary)

产品名称 (Production Name)	TUSC2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUSC2
别名 (Alternative Names)	Tumor suppressor candidate 2, Fusion 1 protein, Fus-1 protein, PDGFA-associated protein 2, TUSC2, C3orf11, FUS1, LGCC, PDAP2
基因 ID (Gene ID)	11334.0
蛋白 ID (SwissProt ID)	O75896. This TUSC2 antibody is generated from a mouse immunized with a recombinant protein of human TUSC2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	12.1kDa

研究背景 (Background)


产品名称: TUSC2 Mouse Monoclonal Antibody
产品货号: AMM86001



May function as a tumor suppressor, inhibiting colony formation, causing G1 arrest and ultimately inducing apoptosis in homozygous 3p21.3 120-kb region-deficient cells.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-TUSC2 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .