

产品名称: CDX1 Mouse Monoclonal Antibody
产品货号: AMM86005



产品概述 (Summary)

产品名称 (Production Name)	CDX1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDX1
别名 (Alternative Names)	Homeobox protein CDX-1, Caudal-type homeobox protein 1, CDX1
基因 ID (Gene ID)	1044.0
蛋白 ID (SwissProt ID)	P47902. Purified His-tagged CDX1 protein was used to produce this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	28.1kDa

研究背景 (Background)


Could play a role in the terminal differentiation of the intestine.

产品名称: CDX1 Mouse Monoclonal Antibody
产品货号: AMM86005



研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-CDX1 Antibody (C-term) at 1:1000-1:2000 dilution

注意事项 (Note)

For research use only .