

产品名称: MKRN2 Mouse Monoclonal Antibody
产品货号: AMM86006



产品概述 (Summary)

产品名称 (Production Name)	MKRN2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MKRN2
别名 (Alternative Names)	Probable E3 ubiquitin-protein ligase makorin-2, 632-, RING finger protein 62, MKRN2, RNF62
基因 ID (Gene ID)	23609.0
蛋白 ID (SwissProt ID)	Q9H000.This MKRN2 antibody is generated from a mouse immunized with recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:500,ICC 1:25-1:50,FC 1:25-1:50
蛋白分子量 (Molecular Weight)	46.9kDa

研究背景 (Background)


产品名称: MKRN2 Mouse Monoclonal Antibody
产品货号: AMM86006



E3 ubiquitin ligase catalyzing the covalent attachment of ubiquitin moieties onto substrate proteins.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-MKRN2 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .