

产品名称: RBM6 Mouse Monoclonal Antibody
产品货号: AMM86007



产品概述 (Summary)

产品名称 (Production Name)	RBM6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBM6
别名 (Alternative Names)	RNA-binding protein 6, Lung cancer antigen NY-LU-12, Protein G16, RNA-binding motif protein 6, RNA-binding protein DEF-3, RBM6, DEF3
基因 ID (Gene ID)	10180.0
蛋白 ID (SwissProt ID)	P78332. This RBM6 antibody is generated from a mouse immunized with a recombinant protein of human RBM6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	128.6kDa

研究背景 (Background)


产品名称: RBM6 Mouse Monoclonal Antibody
产品货号: AMM86007



Specifically binds poly(G) RNA homopolymers in vitro.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-RBM6 Antibody at 1:1000 dilution

注意事项 (Note)

For research use only .