

产品名称: NMI Mouse Monoclonal Antibody
产品货号: AMM86008



产品概述 (Summary)

产品名称 (Production Name)	NMI Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NMI
别名 (Alternative Names)	N-myc-interactor, Nmi, N-myc and STAT interactor, NMI
基因 ID (Gene ID)	9111.0
蛋白 ID (SwissProt ID)	Q13287. This NMI antibody is generated from a mouse immunized with a recombinant protein of human NMI.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:25-1:50, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	35.1kDa

研究背景 (Background)

May be involved in augmenting coactivator protein recruitment to a group of sequence-specific transcription factors.


产品名称: NMI Mouse Monoclonal Antibody
产品货号: AMM86008



Augments cytokine-mediated STAT transcription. Enhances CBP/p300 coactivator protein recruitment to STAT1 and STAT5.

研究领域 (Research Area)

图片 (Image Data)



Anti-NMI Antibody at 1:2000 dilution + A431 whole cell lysate

注意事项 (Note)

For research use only .