

产品名称: S100A2 Mouse Monoclonal Antibody
产品货号: AMM86009



产品概述 (Summary)

产品名称 (Production Name)	S100A2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	S100A2
别名 (Alternative Names)	Protein S100-A2, CAN19, Protein S-100L, S100 calcium-binding protein A2, S100A2, S100L
基因 ID (Gene ID)	6273.0
蛋白 ID (SwissProt ID)	P29034. This S100A2 antibody is generated from a mouse immunized with a recombinant protein of human S100A2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:500, ICC 1:25-1:50
蛋白分子量 (Molecular Weight)	11.1kDa

研究背景 (Background)


产品名称: S100A2 Mouse Monoclonal Antibody
产品货号: AMM86009



May function as calcium sensor and modulator, contributing to cellular calcium signaling. May function by interacting with other proteins, such as TPR-containing proteins, and indirectly play a role in many physiological processes. May also play a role in suppressing tumor cell growth.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-S100A2 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .