

产品名称: RAB18 Mouse Monoclonal Antibody
产品货号: AMM86012



产品概述 (Summary)

产品名称 (Production Name)	RAB18 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB18
别名 (Alternative Names)	Ras-related protein Rab-18, RAB18
基因 ID (Gene ID)	22931.0
蛋白 ID (SwissProt ID)	Q9NP72.This RAB18 antibody is generated from a mouse immunized with a recombinant protein of human RAB18.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:25-1:50
蛋白分子量 (Molecular Weight)	23.0kDa

研究背景 (Background)

Plays a role in apical endocytosis/recycling. May be implicated in transport between the plasma membrane and early


产品名称: RAB18 Mouse Monoclonal Antibody
产品货号: AMM86012



endosomes. Plays a key role in eye and brain development and neurodegeneration.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-RAB18 Antibody at 1:500-1:2000 dilution

注意事项 (Note)

For research use only .