

产品名称: RAB27A Mouse Monoclonal Antibody
产品货号: AMM86013



产品概述 (Summary)

产品名称 (Production Name)	RAB27A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB27A
别名 (Alternative Names)	Ras-related protein Rab-27A, Rab-27, GTP-binding protein Ram, RAB27A, RAB27
基因 ID (Gene ID)	5873.0
蛋白 ID (SwissProt ID)	P51159. This RAB27A antibody is generated from a mouse immunized with a recombinant protein of human RAB27A.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	24.9kDa

研究背景 (Background)


产品名称: RAB27A Mouse Monoclonal Antibody
产品货号: AMM86013



Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

研究领域 (Research Area)

图片 (Image Data)



Anti-RAB27A Antibody at 1:2000 dilution + K562 whole cell lysate

注意事项 (Note)

For research use only .