

产品名称: RRAS2 Mouse Monoclonal Antibody
产品货号: AMM86014



产品概述 (Summary)

产品名称 (Production Name)	RRAS2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RRAS2
别名 (Alternative Names)	Ras-related protein R-Ras2, Ras-like protein TC21, Teratocarcinoma oncogene, RRAS2, TC21
基因 ID (Gene ID)	22800.0
蛋白 ID (SwissProt ID)	P62070. This RRAS2 antibody is generated from a mouse immunized with a recombinant protein of human RRAS2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	24.0kDa

研究背景 (Background)

产品名称: RRAS2 Mouse Monoclonal Antibody
产品货号: AMM86014




It is a plasma membrane-associated GTP-binding protein with GTPase activity. Might transduce growth inhibitory signals across the cell membrane, exerting its effect through an effector shared with the Ras proteins but in an antagonistic fashion.

研究领域 (Research Area)

TGF-beta signaling pathway, MAPK signaling pathway

图片 (Image Data)



All lanes : Anti-RRAS2 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .