

产品名称: RAB5B Mouse Monoclonal Antibody
产品货号: AMM86016



产品概述 (Summary)

产品名称 (Production Name)	RAB5B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB5B
别名 (Alternative Names)	Ras-related protein Rab-5B, RAB5B
基因 ID (Gene ID)	5869.0
蛋白 ID (SwissProt ID)	P61020. This RAB5B antibody is generated from a mouse immunized with a recombinant protein of human RAB5B.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	23.7kDa

研究背景 (Background)


Protein transport. Probably involved in vesicular traffic (By similarity).

产品名称: RAB5B Mouse Monoclonal Antibody
产品货号: AMM86016



研究领域 (Research Area)

图片 (Image Data)



Anti-RAB5B Antibody at 1:2000 dilution + A431 whole cell lysate

注意事项 (Note)

For research use only .