

产品名称: LAP2 Mouse Monoclonal Antibody
产品货号: AMM86017



产品概述 (Summary)

产品名称 (Production Name)	LAP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAP2
别名 (Alternative Names)	Protein LAP2, Densin-180-like protein, Erbb2-interacting protein, Erbin, ERBB2IP, ERBIN, KIAA1225, LAP2
基因 ID (Gene ID)	55914.0
蛋白 ID (SwissProt ID)	Q96RT1. This LAP2 antibody is generated from a mouse immunized with a recombinant protein of human LAP2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	158.3kDa

研究背景 (Background)


产品名称: LAP2 Mouse Monoclonal Antibody
产品货号: AMM86017



Acts as an adapter for the receptor ERBB2, in epithelia. By binding the unphosphorylated 'Tyr-1248' of receptor ERBB2, it may contribute to stabilize this unphosphorylated state.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-LAP2 Antibody at 1:1000-1:2000 dilution

注意事项 (Note)

For research use only .