

产品名称: RASSF2 Mouse Monoclonal Antibody
产品货号: AMM86021



产品概述 (Summary)

产品名称 (Production Name)	RASSF2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF2
别名 (Alternative Names)	Ras association domain-containing protein 2, RASSF2, KIAA0168
基因 ID (Gene ID)	9770.0
蛋白 ID (SwissProt ID)	P50749. This RASSF2 antibody is generated from a mouse immunized with a recombinant protein between 1-326 amino acids from human RASSF2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	37.8kDa

研究背景 (Background)

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes


产品名称: RASSF2 Mouse Monoclonal Antibody
产品货号: AMM86021



STK3/MST2 by protecting it from proteasomal degradation.

研究领域 (Research Area)

图片 (Image Data)



Anti-RASSF2 Antibody at 1:2000 dilution + Raji whole cell lysate

注意事项 (Note)

For research use only .