

产品名称: INA Mouse Monoclonal Antibody
产品货号: AMM86022



产品概述 (Summary)

产品名称 (Production Name)	INA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INA
别名 (Alternative Names)	Alpha-internexin, Alpha-Inx, 66 kDa neurofilament protein, NF-66, Neurofilament-66, Neurofilament 5, INA, NEF5
基因 ID (Gene ID)	9118.0
蛋白 ID (SwissProt ID)	Q16352.This INA antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 290-319 amino acids from human INA.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:200-1:1000,IHC 1:100-1:500,ICC 1:20-1:50
蛋白分子量 (Molecular Weight)	55.4kDa

研究背景 (Background)


产品名称: INA Mouse Monoclonal Antibody
产品货号: AMM86022



Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and they functionally maintain the neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene is a member of the intermediate filament family and is involved in the morphogenesis of neurons.

研究领域 (Research Area)

图片 (Image Data)



Anti-INA Antibody at 1:1000 dilution + SH-SY5Y whole cell lysate

注意事项 (Note)

For research use only .