

产品名称: FAT1 Mouse Monoclonal Antibody
产品货号: AMM86026



产品概述 (Summary)

产品名称 (Production Name)	FAT1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAT1
别名 (Alternative Names)	Protocadherin Fat 1, Cadherin family member 7, Cadherin-related tumor suppressor homolog, Protein fat homolog, Protocadherin Fat 1, nuclear form, FAT1, CDHF7, FAT
基因 ID (Gene ID)	2195.0
蛋白 ID (SwissProt ID)	Q14517. This FAT1 antibody is generated from a mouse immunized with a recombinant protein between 4203-4539 amino acids from human FAT1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	506.3kDa

产品名称: FAT1 Mouse Monoclonal Antibody
产品货号: AMM86026




研究背景 (Background)

Plays an essential role for cellular polarization, directed cell migration and modulating cell-cell contact.

研究领域 (Research Area)

图片 (Image Data)



Anti-FAT1 Antibody at 1:2000 dilution + recombinant protein

注意事项 (Note)

For research use only .