

产品名称: PBRM1 Mouse Monoclonal Antibody
产品货号: AMM86028



产品概述 (Summary)

产品名称 (Production Name)	PBRM1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PBRM1
别名 (Alternative Names)	Protein polybromo-1, hPB1, BRG1-associated factor 180, BAF180, Polybromo-1D, PBRM1, BAF180, PB1
基因 ID (Gene ID)	55193.0
蛋白 ID (SwissProt ID)	Q86U86. This PBRM1 antibody is generated from a mouse immunized with a recombinant protein of human PBRM1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	192.9kDa

研究背景 (Background)


产品名称: PBRM1 Mouse Monoclonal Antibody
产品货号: AMM86028



Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Acts as a negative regulator of cell proliferation.

研究领域 (Research Area)

图片 (Image Data)



Anti-PBRM1 Antibody at 1:1000 dilution + HT-1080 whole cell lysate

注意事项 (Note)

For research use only .