

产品名称: LRP1B Mouse Monoclonal Antibody
产品货号: AMM86030



产品概述 (Summary)

产品名称 (Production Name)	LRP1B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRP1B
别名 (Alternative Names)	Low-density lipoprotein receptor-related protein 1B, LRP-1B, Low-density lipoprotein receptor-related protein-deleted in tumor, LRP-DIT, LRP1B, LRPDIT
基因 ID (Gene ID)	53353.0
蛋白 ID (SwissProt ID)	Q9NZR2. This LRP1B antibody is generated from a mouse immunized with a recombinant protein between 2620-2891 amino acids from human LRP1B.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	515.5kDa

研究背景 (Background)


产品名称: LRP1B Mouse Monoclonal Antibody
产品货号: AMM86030



Potential cell surface proteins that bind and internalize ligands in the process of receptor-mediated endocytosis.

研究领域 (Research Area)

图片 (Image Data)



Anti-LRP1B Antibody at 1:2000 dilution + Recombinant protein

注意事项 (Note)

For research use only .