

产品名称: MLL3 Mouse Monoclonal Antibody
产品货号: AMM86043



产品概述 (Summary)

产品名称 (Production Name)	MLL3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLL3
别名 (Alternative Names)	Histone-lysine N-methyltransferase 2C, Lysine N-methyltransferase 2C, 2.1.1.43, Homologous to ALR protein, Myeloid/lymphoid or mixed-lineage leukemia protein 3, KMT2C, HALR, KIAA1506, MLL3
基因 ID (Gene ID)	58508.0
蛋白 ID (SwissProt ID)	Q8NEZ4. This MLL3 antibody is generated from a mouse immunized with a recombinant protein of human MLL3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	541.4kDa

产品名称: MLL3 Mouse Monoclonal Antibody
产品货号: AMM86043




研究背景 (Background)

Histone methyltransferase. Methylates 'Lys-4' of histone H3. H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation. Central component of the MLL2/3 complex, a coactivator complex of nuclear receptors, involved in transcriptional coactivation. KMT2C/MLL3 may be a catalytic subunit of this complex. May be involved in leukemogenesis and developmental disorder.

研究领域 (Research Area)

图片 (Image Data)



Anti-MLL3 Antibody at 1:2000 dilution + Recombinant protein

注意事项 (Note)

For research use only .