

产品名称: FGFR1OP2 Mouse Monoclonal Antibody
产品货号: AMM86044



产品概述 (Summary)

产品名称 (Production Name)	FGFR1OP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGFR1OP2
别名 (Alternative Names)	FGFR1 oncogene partner 2, FGFR1OP2
基因 ID (Gene ID)	26127.0
蛋白 ID (SwissProt ID)	Q9NVK5. This FGFR1OP2 antibody is generated from a mouse immunized with a recombinant protein between 1-172 amino acids of human FGFR1OP2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	29.4kDa

研究背景 (Background)


May be involved in wound healing pathway.

产品名称: FGFR1OP2 Mouse Monoclonal Antibody
产品货号: AMM86044



研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-FGFR1OP2 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .