

产品名称: ELAVL1 Mouse Monoclonal Antibody
产品货号: AMM86049



产品概述 (Summary)

产品名称 (Production Name)	ELAVL1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ELAVL1
别名 (Alternative Names)	ELAV-like protein 1, Hu-antigen R, HuR, ELAVL1, HUR
基因 ID (Gene ID)	1994.0
蛋白 ID (SwissProt ID)	Q15717. This ELAVL1 monoclonal antibody is generated from mouse immunized with ELAVL1 recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC 1:20-1:50
蛋白分子量 (Molecular Weight)	36.1kDa

研究背景 (Background)

The protein encoded by this gene is a member of the ELAVL


产品名称: ELAVL1 Mouse Monoclonal Antibody
产品货号: AMM86049



protein family. This encoded protein contains 3 RNA-binding domains and binds cis-acting AU-rich elements. It destabilizes mRNAs and thereby regulates gene expression.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti- at 1:1000 dilution

注意事项 (Note)

For research use only .