

产品名称: MYO18B Mouse Monoclonal Antibody
产品货号: AMM86052



产品概述 (Summary)

产品名称 (Production Name)	MYO18B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYO18B
别名 (Alternative Names)	Unconventional myosin-XVIIIb, MYO18B
基因 ID (Gene ID)	84700.0
蛋白 ID (SwissProt ID)	Q8IUG5. This MYO18B antibody is generated from a mouse immunized with recombinant protein from human MYO18B.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:15000, IHC 1:100-1:500, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	285.2kDa

研究背景 (Background)

May be involved in intracellular trafficking of the muscle cell when in the cytoplasm, whereas entering the nucleus, may be


产品名称: MYO18B Mouse Monoclonal Antibody
产品货号: AMM86052



involved in the regulation of muscle specific genes. May play a role in the control of tumor development and progression; restored MYO18B expression in lung cancer cells suppresses anchorage-independent growth.

研究领域 (Research Area)

图片 (Image Data)



Anti-MYO18B Antibody at 1:16000 dilution + Recombinant protein

注意事项 (Note)

For research use only .