

产品名称: RAB20 Mouse Monoclonal Antibody
产品货号: AMM86054



产品概述 (Summary)

产品名称 (Production Name)	RAB20 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB20
别名 (Alternative Names)	Ras-related protein Rab-20, RAB20
基因 ID (Gene ID)	55647.0
蛋白 ID (SwissProt ID)	Q9NX57. This RAB20 antibody is generated from a mouse immunized with recombinant protein from human RAB20.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:4000, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	26.3kDa

研究背景 (Background)

Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf


产品名称: RAB20 Mouse Monoclonal Antibody
产品货号: AMM86054



pathogens, such as *S.aureus* and *M.tuberculosis*. Plays a role in the fusion of phagosomes with lysosomes.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-RAB20 Antibody at 1:2000-1:4000 dilution

注意事项 (Note)

For research use only .