

产品名称: IFNB1 Mouse Monoclonal Antibody
产品货号: AMM86055



产品概述 (Summary)

产品名称 (Production Name)	IFNB1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFNB1
别名 (Alternative Names)	Interferon beta, IFN-beta, Fibroblast interferon, IFNB1, IFB, IFNB
基因 ID (Gene ID)	3456.0
蛋白 ID (SwissProt ID)	P01574. This IFNB1 antibody is generated from a mouse immunized with a recombinant protein between 1-187 amino acids from human IFNB1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	22.3kDa

研究背景 (Background)

Has antiviral, antibacterial and anticancer activities.


产品名称: IFNB1 Mouse Monoclonal Antibody
产品货号: AMM86055



研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Anti-IFNB1 Antibody at 1:2000 dilution + IFNB1 protein lysate

注意事项 (Note)

For research use only .