

产品名称: MCFD2 Mouse Monoclonal Antibody
产品货号: AMM86057



产品概述 (Summary)

产品名称 (Production Name)	MCFD2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCFD2
别名 (Alternative Names)	Multiple coagulation factor deficiency protein 2, Neural stem cell-derived neuronal survival protein, MCFD2, SDNSF
基因 ID (Gene ID)	90411.0
蛋白 ID (SwissProt ID)	Q8NI22.This MCFD2 monoclonal antibody is generated from mouse immunized with MCFD2 recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	16.4kDa

研究背景 (Background)


产品名称: MCFD2 Mouse Monoclonal Antibody
产品货号: AMM86057



The MCFD2-LMAN1 complex forms a specific cargo receptor for the ER-to-Golgi transport of selected proteins. Plays a role in the secretion of coagulation factors.

研究领域 (Research Area)

图片 (Image Data)



Anti-MCFD2 Antibody at 1:1000 dilution + 293T/17 whole cell lysate

注意事项 (Note)

For research use only .