

产品名称: THEM2 Mouse Monoclonal Antibody
产品货号: AMM86059



产品概述 (Summary)

产品名称 (Production Name)	THEM2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THEM2
别名 (Alternative Names)	Acyl-coenzyme A thioesterase 13, Acyl-CoA thioesterase 13, 312-, Thioesterase superfamily member 2, Acyl-coenzyme A thioesterase 13, N-terminally processed, ACOT13, THEM2
基因 ID (Gene ID)	55856.0
蛋白 ID (SwissProt ID)	Q9NPJ3.This THEM2 monoclonal antibody is generated from mouse immunized with THEM2 recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	15.0kDa

产品名称: THEM2 Mouse Monoclonal Antibody
产品货号: AMM86059




研究背景 (Background)

THEM2 belongs to the UPF0152 (paal) family. Its exact function is still undetermined.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-THEM2 Antibody at 1:1000 dilution

注意事项 (Note)

For research use only .