

产品名称: RAB21 Mouse Monoclonal Antibody
产品货号: AMM86061



产品概述 (Summary)

产品名称 (Production Name)	RAB21 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat, African Green Mouseonkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB21
别名 (Alternative Names)	Ras-related protein Rab-21, RAB21, KIAA0118
基因 ID (Gene ID)	23011.0
蛋白 ID (SwissProt ID)	Q9UL25. This RAB21 antibody is generated from a mouse immunized with a recombinant protein between 1-225 amino acids from the human RAB21.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	24.3kDa

研究背景 (Background)

Regulates integrin internalization and recycling, but does not influence the traffic of endosomally translocated receptors in


产品名称: RAB21 Mouse Monoclonal Antibody
产品货号: AMM86061



general. As a result, may regulate cell adhesion and migration (By similarity). During the mitosis of adherent cells, controls the endosomal trafficking of integrins which is required for the successful completion of cytokinesis. Involved in neurite growth (By similarity).

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-RAB21 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .