

产品名称: CBFA2T3 Mouse Monoclonal Antibody
产品货号: AMM86062



产品概述 (Summary)

产品名称 (Production Name)	CBFA2T3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CBFA2T3
别名 (Alternative Names)	Protein CBFA2T3, MTG8-related protein 2, Myeloid translocation gene on chromosome 16 protein, hMTG16, Zinc finger MYND domain-containing protein 4, CBFA2T3, MTG16, MTGR2, ZMYND4
基因 ID (Gene ID)	863.0
蛋白 ID (SwissProt ID)	O75081. This CBFA2T3 antibody is generated from a mouse immunized with a recombinant protein between 418-653 amino acids from human CBFA2T3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	71.2kDa

产品名称: CBFA2T3 Mouse Monoclonal Antibody
产品货号: AMM86062




研究背景 (Background)

Functions as a transcriptional repressor. Regulates the proliferation and the differentiation of erythroid progenitors by repressing the expression of TAL1 target genes. Plays a role in granulocyte differentiation.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-CBFA2T3 Antibody at 1:500-1:2000 dilution

注意事项 (Note)

For research use only .