

产品名称: GOLGA5 Mouse Monoclonal Antibody
产品货号: AMM86063



产品概述 (Summary)

产品名称 (Production Name)	GOLGA5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GOLGA5
别名 (Alternative Names)	Golgin subfamily A member 5, Cell proliferation-inducing gene 31 protein, Golgin-84, Protein Ret-II, RET-fused gene 5 protein, GOLGA5, RETII, RFG5
基因 ID (Gene ID)	9950.0
蛋白 ID (SwissProt ID)	Q8TBA6. This GOLGA5 antibody is generated from a mouse immunized with a recombinant protein between 1-300 amino acids from human GOLGA5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	83.0kDa

研究背景 (Background)


产品名称: GOLGA5 Mouse Monoclonal Antibody
产品货号: AMM86063



Involved in maintaining Golgi structure. Stimulates the formation of Golgi stacks and ribbons. Involved in intra-Golgi retrograde transport.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-GOLGA5 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .