

产品名称: FGFR1OP2 Mouse Monoclonal Antibody
产品货号: AMM86065



产品概述 (Summary)

产品名称 (Production Name)	FGFR1OP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGFR1OP2
别名 (Alternative Names)	FGFR1 oncogene partner 2, FGFR1OP2
基因 ID (Gene ID)	26127.0
蛋白 ID (SwissProt ID)	Q9NVK5. This FGFR1OP2 antibody is generated from a mouse immunized with a recombinant protein between 1-172 amino acids of human FGFR1OP2.

产品应用 (Application)

稀释比 (Dilution Ratio)	FC 1:25-1:50
蛋白分子量 (Molecular Weight)	29.4kDa

研究背景 (Background)


May be involved in wound healing pathway.

产品名称: FGFR1OP2 Mouse Monoclonal Antibody
产品货号: AMM86065



研究领域 (Research Area)

图片 (Image Data)



Overlay histogram showing U-2OS cells stained with AM8701b(green line). The cells were fixed with 2% paraformaldehyde (10 min) and then permeabilized with 90% methanol for 10 min. The cells were then incubated in 2% bovine serum albumin to block non-specific protein-protein interactions followed by the antibody (AM8701b, 1:25 dilution) for 60 min at 37°C. The secondary antibody used was Goat-Anti-Mouse IgG, DyLight® 488 Conjugated Highly Cross-Adsorbed(OJ192088) at 1/200 dilution for 40 min at 37°C. Isotype control antibody (blue line) was mouse IgG2a (1µg/1x10⁶ cells) used under the same conditions. Acquisition of >10, 000 events was performed.

注意事项 (Note)

For research use only .