

产品名称: PFKP Mouse Monoclonal Antibody
产品货号: AMM86066



产品概述 (Summary)

产品名称 (Production Name)	PFKP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, African Green Mouseonkey, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PFKP ATP-dependent 6-phosphofructokinase, platelet type {ECO:0000255 HAMAP-Rule:MF_03184}, ATP-PFK {ECO:0000255 HAMAP-Rule:MF_03184}, PFK-P, 2.7.1.11 {ECO:0000255 HAMAP-Rule:MF_03184}, 6-phosphofructokinase type C, Phosphofructo-1-kinase isozyme C, PFK-C, Phosphohexokinase {ECO:0000255 HAMAP-Rule:MF_03184}, PFKP, PFKF
别名 (Alternative Names)	
基因 ID (Gene ID)	5214.0
蛋白 ID (SwissProt ID)	Q01813.This PFKP antibody is generated from a mouse immunized with a recombinant protein between 40-300 amino acids from the C-terminal region of human PFKP.

产品应用 (Application)

产品名称: PFKP Mouse Monoclonal Antibody
产品货号: AMM86066




稀释比 (Dilution Ratio) WB 1:1000-1:2000
蛋白分子量 (Molecular Weight) 85.6kDa

研究背景 (Background)

Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of glycolysis.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-PFKP Antibody (C-term) at 1:4000 dilution

注意事项 (Note)

For research use only .