

产品名称: PACSIN2 Mouse Monoclonal Antibody
产品货号: AMM86072



产品概述 (Summary)

产品名称 (Production Name)	PACSIN2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PACSIN2
别名 (Alternative Names)	Protein kinase C and casein kinase substrate in neurons protein 2, Syndapin-2, Syndapin-II, PACSIN2
基因 ID (Gene ID)	11252.0
蛋白 ID (SwissProt ID)	Q9UNF0. This PACSIN2 antibody is generated from a mouse immunized with a recombinant protein between 250-486 amino acids from human PACSIN2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	55.7kDa

研究背景 (Background)


产品名称: PACSIN2 Mouse Monoclonal Antibody
产品货号: AMM86072



Lipid-binding protein that is able to promote the tubulation of the phosphatidic acid-containing membranes it preferentially binds. Plays a role in intracellular vesicle- mediated transport. Involved in the endocytosis of cell-surface receptors like the EGF receptor, contributing to its internalization in the absence of EGF stimulus. May also play a role in the formation of caveolae at the cell membrane. Recruits DNM2 to caveolae, and thereby plays a role in caveola-mediated endocytosis.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-PACSIN2 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .