

产品名称: NPPA Mouse Monoclonal Antibody
产品货号: AMM86073



产品概述 (Summary)

产品名称 (Production Name)	NPPA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPPA
别名 (Alternative Names)	Natriuretic peptides A, CDD-ANF, Cardiodilatin, CDD, Cardiodilatin-related peptide, CDP, Prepronatriodilatin, Atrial natriuretic factor, ANF, Atrial natriuretic peptide, ANP, NPPA, ANP, PND
基因 ID (Gene ID)	4878.0
蛋白 ID (SwissProt ID)	P01160. This NPPA antibody is generated from a mouse immunized with a recombinant protein of human NPPA.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	16.4kDa

产品名称: NPPA Mouse Monoclonal Antibody
产品货号: AMM86073




研究背景 (Background)

Hormone playing a key role in cardiovascular homeostasis through regulation of natriuresis, diuresis, and vasodilation. Also plays a role in female pregnancy by promoting trophoblast invasion and spiral artery remodeling in uterus. Specifically binds and stimulates the cGMP production of the NPR1 receptor. Binds the clearance receptor NPR3.

研究领域 (Research Area)

图片 (Image Data)



Anti-NPPA Antibody at 1:2000 dilution + mouse heart lysate

注意事项 (Note)

For research use only .