

产品名称: DCLK1 Mouse Monoclonal Antibody
产品货号: AMM86074



产品概述 (Summary)

产品名称 (Production Name)	DCLK1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DCLK1
别名 (Alternative Names)	Serine/threonine-protein kinase DCLK1, 2.7.11.1, Doublecortin domain-containing protein 3A, Doublecortin-like and CAM kinase-like 1, Doublecortin-like kinase 1, DCLK1, DCAMKL1, DCDC3A, KIAA0369
基因 ID (Gene ID)	9201.0
蛋白 ID (SwissProt ID)	O15075. This DCLK1 antibody is generated from a mouse immunized with a recombinant protein of human DCLK1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	82.2kDa

产品名称: DCLK1 Mouse Monoclonal Antibody
产品货号: AMM86074




研究背景 (Background)

Probable kinase that may be involved in a calcium- signaling pathway controlling neuronal migration in the developing brain. May also participate in functions of the mature nervous system.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-DCLK1 Antibody at 1:4000 dilution

注意事项 (Note)

For research use only .